





A NEW REVOLUTION IN CELL COUNTING STARTS HERE

iWash® Slide Cleaner can now clean most of the cell counter slides!



NOW IT IS POSSIBLE TO WASH AND RE-USE CELL COUNTER SLIDES IN LESS THAN 20 SECONDS!

- Less than 15 seconds to wash each slide
- Reuse each slide >20 times
- Reduce the cost of cell counting >95%
- Turn your cell counter into almost Slide-free operation
- Reduce the single-use plastic waste and become green
- Worldwide proven technology with fantastic testimonials
- Compatible with all major automated cell counters in the World

THE IWASH® SLIDE CLEANER FROM IMRALI INVENTIONS A NEW REVOLUTION IN CELL COUNTING

JOIN THE REVOLUTION!

SAVE MONEY, REDUCE WASTE, INCREASE PRODUCTIVITY AND IMPROVE SUSTAINABILITY.

Specially designed Accesories:




HumiBox ®: Stops slides from drying out upto 48 hours, if you prefer to wash them later.



Slide holding rack: Holds the washed slides organised in one place, keeps them scratch and dirt free.

New versus washed slide comparision



The iWash® Slide Cleaner is perfect for the laboratory that wants to minimise the cost of automated cell counting and to reduce hazardous waste. Reducing costs and improving sustainability are major issues within most facilities and the iWash® Slide Cleaner was designed after years of development in partnership with the scientific community to meet these needs.

iWash® is easily integrated into your cell counting workflow and with our patented direct injection technology it can quickly and efficiently cleans slides without compromise so it is ready for the next analysis in less than 15 seconds.

WELL PROVEN TECHNOLOGY, We have got fantastic testimonials from all around the world from USA, Brazil, EU, China, Singapore, Cyprus, Turkey, Pakistan, Australia and many more...

For more information or to arrange a demo please contact us.

We are the inventor and manufacturer of the iWash® systems

E mail: info@imraliinventions.co.uk
Website: www.imrali.co.uk





Used all around the world